

seT



(seT capsules)

SELF ETCHING, SELF ADHESIVE RESIN CEMENT

SELBSTÄTZENDER, SELBSTKLEBENDER KUNSTSTOFF-ZEMENT

CIMENTO RESINOSO AUTO-CONDICIONANTE E AUTO-ADESIVO

CEMENTO DE RESINA, AUTOGRABANTE Y AUTOADHESIVO

RÉSINE DE SCELLEMENT AUTO-ADHÉSIVE AUTO-MORDANÇANT

CEMENTO COMPOSITO AUTOMORDENZANTE ED AUTOADESIVO

ZELF ETSEND, ZELF ADHESIEF RESIN CEMENT

SELVÆTSENDE, SELVADHÆRERENDE RESIN CEMENT

SELVETSENDE, SELVBINDENDE RESIN SEMENT

ITSE-ETSAAVA, ITSEKIINNITTYVÄ MUOVISEMENTTI

ΑΥΤΟΑΔΡΟΠΟΙΟΥΜΕΝΗ, ΑΥΤΟΣΥΓΚΟΛΛΟΥΜΕΝΗ ΡΗΤΙΝΩΔΗΣ ΚΟΝΙΑ

SAMOWYTRAWIAJĄCY CEMENT KOMPOZYTOWY NIE WYMAGAJĄCY UŻYCIA SYSTEMU WIĄŻĄCEGO

ÖNSAVAZÓ ÉS ÖNKÖTŐ MŰGYANTABÁZISÚ RAGASZTÓCEMENT

ISESÖÖVITAV, ISEADHESEERUV TSEMENT

SAMOJEDKAJOČ, SAMOADHEZIVEN SMOLNAT CEMENT

PAŠKODINOŠS, PAŠADHEZĪVS SVEĶU CEMENTS

SAVAIME ĒSDINANTIS BEI SURIŠANTIS DERVŪ CEMENTAS

SAMOLEPTACÍ, ADHEZIVNÍ CEMENT ZESÍLENÝ PRYSKYŘICÍ

SAMOLEPTACÍ, ADHEZÍVNÝ, ŽIVICOVÝ CEMENT

歯科接着用レジンセメント セルフエッチングデュアルキュアタイプ

自酸蚀、自粘结树脂粘剂





SDI's top rated cement, seT now available in syringes


seT



seT is a self etching, self adhesive, resin cement, designed for final cementation of metal, resin, all ceramic and ceramic based restoratives. seT chemically bonds to tooth structure and all types of core material. seT is available in capsules and a paste/paste automix syringe system (seT PP).

Proprietary ionglass™ filler
seT utilises SDI's proprietary *ionglass™* filler developed by our glass technologists. *ionglass™* is a radiopaque, high ion releasing, reactive glass used in SDI's range of dental cements. seT releases substantially higher fluoride to assist with remineralisation of the natural dentition.




High fluoride release
seT's free movement of fluoride provides benefits to the tooth. Fluoride plays several significant roles in any caries - prevention program. These include the formation of fluorapatite, which is more acid resistant than hydroxyapatite. Fluorapatite is harder than other tooth minerals. It helps to provide protection from tooth decay.



Fast and easy

seT is an all in one luting cement. It etches, bonds and fills with just the one component. Simply apply seT and seat your restoration.

seT produces a low film thickness and fluidity to seat cemented restorations.



No Bisphenol A or HEMA

seT does not contain any Bisphenol A (including its derivatives) or HEMA. Use this product on your patients with confidence and peace of mind.

Dual cure

Designed to combine the advantages of light curing and self curing, seT is the ideal luting cement for most of your cementation needs[^].

Radiopacity

Detectable on x-ray radiographs, seT is great for recalls when margins are regularly checked for integrity.

Minimal to no post-operative sensitivity

seT does not require separate etching and washing. The smear layer remains intact to protect the tooth. This dramatically reduces, if not eliminates, post-operative sensitivity.

Choice of delivery tips

seT PP syringes can be adapted for a variety of delivery options. The regular mixing tips are for predominantly used for dispensing into crown and bridge. In conjunction with a wide mixing tip, endodontic and intraoral tips are also available.



Extensive colour range

seT is available in 5 shades: Translucent, A1, A2, OA3 and Opaque White. The large variety of shades helps maximise aesthetics.

Indications

- Cementation of metal, resin, all ceramic and ceramic based inlays, onlays[^], bridges and crowns.
- Cementation pins, posts and screws

[^] In the case of veneers, use an adhesive (Go!, Stae, etc according to manufacturer's instructions including any etching procedure) prior to seT to improve the retentive strength. Not suitable for veneers requiring extended working time.

**seT PP.

instructions:

- 1 Remove temporary cement residue and rinse thoroughly. Remove excess water, keep moist.



- 2 Prepare restoration according to manufacturer's instructions, and check fit.



- 3 Extrude cement evenly onto all bonding surfaces of the restoration and tooth.



- 4 Seat restoration, and leave for 30 seconds. Tack cure excess for 2 seconds or wait 2 minutes.



- 5 Trim excess material away.



- 6 Maintain isolation until set is verified (approximately 5 minutes), or light cure through translucent restoration for at least 20 seconds.



seT



seT PP automix syringe refill

2 seT PP 7g automix syringes
28 automix syringe mixing tips

Reorder:

8805101 A1
8805102 A2
8805103 OA3
8805104 White Opaque
8805105 Translucent

seT PP assorted automix syringe kit

1 each of A1, A2, OA3, White Opaque and Translucent
7g seT PP automix syringes
70 automix syringe mixing tips
Reorder 8810002

seT capsules assorted kit

10 each of A1, A2, OA3, White Opaque & Translucent
seT capsules
Reorder: 8810000

seT capsules

50 seT capsules
Reorder:
8805001 A1
8805002 A2
8805003 OA3
8805004 White Opaque
8805005 Translucent

50 x seT PP Mixing Tips - Regular

Reorder: 8810003

10 x seT PP Mixing Tips - Wide

Reorder: 8810004

10 x seT PP Oral Tips

(To be used in conjunction with seT PP Mixing Tips - Wide)
Reorder: 8810005

10 x seT PP Endo Tips

(To be used in conjunction with seT PP Mixing Tips - Wide)
Reorder: 8810006



Riva Applicator 2[†]
Reorder 5545013

Riva Applicator[†]
Reorder 5545009

*Not the registered trademarks of SDI limited.

(1) Source: Al-Naimi O.T., McCabe J. F. Fluoride release of seT self-etch/self-adhesive resin cement with that of Competitor Products. (18 months report) Newcastle University, U.K. February 2007.

(2) SDI Test Data.

seT: THE DENTAL ADVISOR, Vol. 25, No. 08. October 2008

seT PP: THE DENTAL ADVISOR, Vol. 27, No. 7, Sept. 2010, THE DENTAL ADVISOR preferred product 2011.



Made in Australia by SDI Limited
Bayswater, Victoria 3153
Australia 1 800 337 003
Austria 00800 022 55 734
Brazil 0800 770 1735
France 00800 022 55 734
Germany 0800 100 5759
Ireland 01 886 9577
Italy 800 780625
New Zealand 0800 734 034
Spain 00800 022 55 734
United Kingdom 00800 022 55 734
USA & Canada 1 800 228 5166
www.sdi.com.au



M880551 E

+0036188055118