

# **ACKURETTA UV BOX**



The Ackuretta UV Box is a compact curing box lined with long-life ultraviolet LED bulbs for enhanced curing performance. The reflective stainless steel chamber and bottom mirror surface provides even curing to all sides. The UV Oven is designed with two modes, continuous and flash (10 flashes/second), to ensure all resins that react to wavelengths 365, 385, and 405nm can cure completely and efficiently, with strengthened print hardness. It is also suitable for Class I materials.



**Small carbon footprint** 



**Mirror reflection** enhances curing performance



365/385/405 nm mix matrix LEDs for a thorough cure



**Flash function** 10 flashes/sec strengthens print hardness



Suitable for Class I materials

## **UV BOX PACKAGE INCLUDES**

Main curing box equipment 1 Pcs
1 Power cord 1 Pcs
1 Protective Goggles 1 Pair
1 User Manual 1 Pcs

# **TECHNICAL SPECIFICATIONS**



Power Supply Voltage 100-240 VAC 50/60 Hz



Wattage Max 65 W



Dimensions H 150 x W 200 x D 216mm



Weight 3.3Kg

### **Operating Environment**



Ambient Temperature  $+10^{\circ}\text{C} \sim +40^{\circ}\text{C}$ 



Relative Humidity  $0\% \sim 75\%$ 



Ambient Pressure 700hPa  $\sim$  1060hPa

#### Transport & Storage Condition



Ambient Temperature  $-10^{\circ}\text{C} \sim +70^{\circ}\text{C}$ 



Relative Humidity  $10\% \sim 90\%$ 



Ambient Pressure 500hPa ∼ 1060hPa

